

Thermal path optimization between cryocooler and sample

Marcel van Wensveen

Rens van Alebeek


Ron van den Boogaard

Ronald Lamers

13 March 2024

euspen conference

 The world leader in serving science



Contents

- Introduction
- Design and components
- Thermal modelling
- Verification measurements
- Improvements
- Results and conclusions


Introduction

- Thermo Fisher Scientific Eindhoven develops and produces Transmission Electron Microscopes (TEM).
- In a new development, TEM and Atom Probe Tomography (APT) will be combined.
- Existing TEM requirements:
 - Tight volume (in the heart of the TEM)
 - Only non-magnetic materials allowed
- For APT, some new system requirements are added:
 - Sample temperature < 35 Kelvin
 - Ultra-High Vacuum around the sample (1E^{-10} mbar)
 - The sample needs to be at high voltage (12 kV)
- A thermal connection between cooler and sample needed to be realized.




Design and main components

- The main components of the design:




Design and main components

- Thermal requirements:
 - Temperature delta
 - Heat load



Design – Emissivity & Coatings

- Heat load needs to be minimized.
- Low emissivity (gold) coating.
- Shield with high emissivity coating.



Modelling – Non-linearities


- FE modelling has been used to optimize the design.

- At cryogenic temperatures:

- Limited data sources
- Material properties highly non-linear


- Residual-resistance ratio of copper

- Purity
- Crystallographic defects
- Heat treatment
- Cold deformation
- Dimensions




Modelling – Results

- After several design iterations, a design within specification has been obtained.
- The main temperature gradients are over the braid and stage bulk material.





Measurements - Cryostat

- To verify the model and design, tests have been performed.
- A cryostat has been used to perform cryogenic tests.




Measurements – Thermal conductivity

- For multiple material, the thermal conductivity has been measured.
- Below, a low thermal conductivity titanium alloy test setup is shown.




Measurement – Results

- All components have been tested.
- Temperature gradients due to material conductivity are as expected.
- Temperature gradients over contacts are much higher.



Thermal contact resistance

- Equations of Mikic* have been used.
- Some resistances where a factor ~ 10.000 off.
- At cryogenic temperature, no realiable equations.
- Many measurements have been performed.
- Improvement options are available:
 - Increase contact force
 - Applying an intermediate layer or coating
 - Optimizing surface roughness & flatness
 - Soldering
 - Welding




Source: Comsol

*B.B. Mikic, Thermal contact conductance; theoretical considerations, 1974, International Journal of Heat and Mass Transfer. 17.2: 205-217


Thermal contact resistance

- Multiple test setups for thermal contact resistance.
- Test performed with:
 - Contact force variations
 - (Thick) gold coating
 - Tin and indium interface layers
 - Soldering / welding



Final results

- After all improvements of the contacts have been applied, the following results have been obtained.
- Soldering and indium gave the main improvements



Summary and conclusions

- Many thermal challenges
 - Material properties
 - Thermal contact resistance
 - Testing is of major importance
- System within the specification!



End

Thank you for listening!

